



ORTHOPEDIC & SPINE INSTITUTE

at St. Mary's Medical Center

Orthex Computer Dependent External Fixation & Drive Rail Monolateral External Fixation

Dror Paley, MD, FRCSC

CEO, Founder & Director

Paley Orthopedic & Spine Institute

St. Mary's Medical Center, West Palm Beach, Florida, USA

Professor of Orthopedics, Florida Atlantic University



ORTHEX Frame


Ring Design

Two rows:


- inner row holes
- outer row slots

Advantages:


- *More Bending/Cantilever strength*
- *More options for pin/wire/strut fixation*



Concentric Ring Design



- Inner row matches next size down outer row
- Outer row matches next size up inner row outer row slots



Advantages:

- *Can connect either via threaded rods or struts*
- *Tapered frames to fit anatomy when using struts*
- *Cylindrical frames when using threaded rods*


Special Rings

- Foot Rings
- 2/3 Rings
- Low Profile Arch



Closing Foot Ring




- Foot Rings: short and long
- 2/3 Rings
- Low Profile Arch



Adult & Pediatric Rings

Full Rings, Partial Ring, Footplate, Specialized Arch



- **Adult:** 105, 130, 155, 180, 205, and 230mm*
- **Small-Bone/Pediatric:** XS, Small, Medium, Large
- **Peds frame is 60% lighter**



Adult & Pediatric Struts


Full Rings, Partial Ring, Footplate, Specialized Arch

- **Adult:** 105, 130, 155, 180, 205, and 230mm*
- **Small-Bone/Pediatric:** XS, Small, Medium, Large
- **Peds frame is 60% lighter**



Strut Design

Ball and Socket Construction



Paley Frame: No Vibration

TSF : Vibration



TSF-vibration



Orthex-no vibration

Struts

- Telescoping Struts
- Improved Range
- Fewer Sizes,
- Fewer Strut Changes



Fully collapsed strut, 137mm



Partially Extended, 215mm




Fully Extended, 295mm

3 Strut Sizes



All Struts are both Standard & Fast Fix



TSF

	Standard Struts				Fast Fx ©			
Size:	Extra Short	Short	Medium	Long	Extra Short	Short	Medium	Long
Range:	75-96 mm	90-125 mm	116-178 mm	170-280 mm	91-121mm	116-152 mm	143-205 mm	195-311 mm
Cat. No.:	7107-0205	7107-0210	7107-0220	7107-0230	7107-0705	7107-0710	7107-0720	7107-0730

Orthex struts are the most efficient

Efficiency: max-min / min

Orthex PaleyFrame				
Size	Min.	Max	Range	Efficiency
Small	70	105	35	0.5
Medium	93	170	77	0.8
Large	137	295	158	1.2

*Z-Plate Strut Extender			
Z-Min	Z-Max	Range	Efficiency
20	155	135	6.8
43	220	177	4.1
87	345	258	3.0

Smith & Nephew TSF				
Size	Min.	Max	Range	Efficiency
Extra Small	91	121	30	0.3
Small	116	152	36	0.3
Medium	143	205	62	0.4
Large	195	311	116	0.6

Orthofix TL-Hex				
Size	Min.	Max.	Range	Efficiency
Small	92	122	30	0.3
Medium	114	184	70	0.6
Large	158	318	160	1.0

Strut with Z-Plate



Fewer Strut Changes

- Fewer Strut Changes
- Reduced Inventory



Strut Appointments			
Orthex	Adjust Only	Change Only	Smith & Nephew
4	2	2	4
4	0	4	6
2	0	2	5
7	4	3	7
6	4	2	5
0	0	0	3
4	3	1	5
4	2	2	6
6	1	5	9
0	0	0	3
37	16	21	53

Z-plates allow Extreme frame angles: 90° angle with rings touching





Z-plates: extend the longest strut and shorten the shortest strut




Z-plate extension



Extreme Frame



Extreme Frame




During Correction




Half-Pin Blocks

- Double column half pin blocks
- Pin clamp effect effect due to wider pin grab (strengthens the cantilever of the pin)
- More stable/prevents pins bending/need fewer pins
- Improved drill guidance and accuracy





Half-Pin Blocks

- One hole to five hole peds and adult




PaleyPin™

- Variable Pitch Thread Half Pin (pat pend)
- Extended HA Coated Half Pin




- Hole Locator Tip
- Tapered Core/ Non-tapered Threads
 - Adult (3.0mm, 4.5mm, 6.0mm OD) (10-60mm)
 - Peds (3.0mm, 4.0mm, 5.0mm OD) (10-50)

HA Coated Pins & Wires






PaleyPin™




Ring Positioning Software (RPS)




Fibular Hemimelia




50mm Lengthening 20 deg Procurvatum




CFD & FH



Foot Frame





Foot Frame





Distal Radial Growth Arrest




Radial Growth Arrest





Pediatric Frame - 18mo. old child



Severe Genu Valgum & Procurvatum Bilateral (Fanconi Syndrome)



Two Level Acute Correction Femurs Gradual Correction Tibias




10-29-2018

Two Level Acute Correction Femurs Gradual Correction Tibias



11-12-2018




11-16-2018

Full Correction Achieved





2-4-2019

After Lengthening Completed: removal



3 /4/2019



3/5/2019

Final Result: Fracture Left Tibia with Varus




8/5/2019

Drive Rail

A New Level in Rail Technology



Modular Rail Segments



Rail Assembly




Unique Lengthening Mechanism (patented)

no distraction rod



Assembly of Lengthening Mechanism



Assembly of Cars (pin clamps)




End Car




Wide Variety of Cars: Pin Clamps





Lengthening




Basic Assembly



Add a Hinge



Hinge Options on End Car/Rail




Hinge option on Moving Car e.g. Hip Hinge




Extend Rail from Hinge





Add an Arch




Arch Connectors




Hip and Knee Hinge Construct




Hip & Knee Hinges with Lengthening




Angular Correction with Swivel Car





Swivel Car Connected to Knee Hinge



First Clinical Use of Drive Rail July 2015




Lengthening Femur & Tibia



Lengthening Arms: Dwarfism




Lengthening Femur with Knee Hinge




CFD Lengthening




Elbow Contracture



CFD: No Hip; Very Short Femur Ex Fix: Pelvis-Femur-Tibia









CFD Lengthening





CFD Lengthening




Tibial Hemimelia



Forearm Lengthening





Summary

- Most Modular Rail System
- Can *Span Joints* better than any other Rail
- Peds and Adult with one Rail size
- Unique Lengthening mechanism
- Most efficient design
- Integrates well with Circular Device